

Package ‘MiChip’

August 21, 2025

Version 1.62.0

Date 2009-09-10

Title MiChip Parsing and Summarizing Functions

Author Jonathon Blake <blake@embl.de>

Maintainer Jonathon Blake <blake@embl.de>

Depends R (>= 2.3.0), Biobase

Imports Biobase

Description This package takes the MiChip miRNA microarray .grp scanner output files and parses these out, providing summary and plotting functions to analyse MiChip hybridizations. A set of hybridizations is packaged into an ExpressionSet allowing it to be used by other BioConductor packages.

License GPL (>= 2)

biocViews Microarray, Preprocessing

git_url <https://git.bioconductor.org/packages/MiChip>

git_branch RELEASE_3_21

git_last_commit f99b748

git_last_commit_date 2025-04-15

Repository Bioconductor 3.21

Date/Publication 2025-08-20

Contents

boxplotData	2
boxplotDataNoFile	3
correctForFlags	3
MiChip	4
myForgivingMedian	4
naOmitMedian	5
normalizePerChipMedian	5
outputAnnotatedDataMatrix	6

panelCor	6
parseRawData	7
plotIntensitiesScatter	8
removeUnwantedRows	8
returnAnnotatedDataMatrix	9
setIntensityCutoff	9
standardRemoveRows	10
summarizeIntensitiesAsMedian	11
workedExampleMedianNormalize	11
workedExampleNotNormalizedData	12
Index	13

boxplotData	<i>Create Boxplot of data</i>
-------------	-------------------------------

Description

Creates a boxplot of expression data contained in a matrix and writes this to a file.

Usage

boxplotData(dmat, exptname, dlevel)

Arguments

- | | |
|----------|--|
| dmat | matrix containing expression data to be boxplotted |
| exptname | Name of the experiment, used to build filename |
| dlevel | Stage of the experiment e.g. raw, summarized, normalized |

Examples

```
#Create a file of a boxplot containing normalized expression data for myexpt
## Not run:
boxplotData(dmat, "MyExpt", "mednormed")

## End(Not run)
```

boxplotDataNoFile	<i>Create Boxplot of data</i>
-------------------	-------------------------------

Description

Creates a boxplot of expression data contained in a matrix.

Usage

```
boxplotDataNoFile(dmat, exptname, dlevel)
```

Arguments

dmat	matrix containing expression data to be boxplotted
exptname	Name of the experiment, used to build plot title
dlevel	Stage of the experiment e.g. raw, summarized, normalized

Examples

```
#Create a boxplot containing normalized expression data for myexpt
## Not run: boxplotDataNoFile(dmat, "MyExpt", "mednormed")
```

correctForFlags	<i>Corrects for spots flagged as not present</i>
-----------------	--

Description

Spots flagged with a -ve quality flag value by the scanner may be regarded as not present. This method sets their intensity to NA.

Usage

```
correctForFlags(eset, intensityCutoff=0)
```

Arguments

eset	ExpressionSet containing intensity values and flags to be filtered
intensityCutoff	value of lowest acceptable intensity value in the experiment

Examples

```
#Correct ExpressionSet for flags of spots marked as unreadable
## Not run:
myCorrectedEset <-correctForFlags(eset, intensityCutoff=0)

## End(Not run)
```

MiChip

Introduction to the MiChip Package

Description

A library for processing MiChip hybridizations

Author(s)

Jonathon Blake

myForgivingMedian

Produce Median from Probe Intensity values

Description

Creates a median to summarize the intensities for individual probes, giving that not all probes will have a valid intensity reading

Usage

```
myForgivingMedian(mat, minSumlength=0)
```

Arguments

mat	matrix of data to calculate the median from
minSumlength	The lowest acceptable length of the matrix to calculate a median

Examples

```
#Calculate the median of a matrix omitting NAs
## Not run:
myForgivingMedian(mat, minSumlength=0)

## End(Not run)
```

naOmitMedian	<i>Calculates the median of an array excluding NAs</i>
--------------	--

Description

Calculates the median of an array excluding NAs

Usage

```
naOmitMedian(mat, madAdjust=FALSE)
```

Arguments

mat	A single dimensional matrix
madAdjust	if TRUE then summarized data will be filtered according to the MAD median absolute deviation and set to NA if the median is less than MAD

Examples

```
#Calculate the median of matrix mat omiting NAs
## Not run:
myMedian <-naOmitMedian(mat, madAdjust=TRUE)

## End(Not run)
```

normalizePerChipMedian	<i>Normalize to median intensity</i>
------------------------	--------------------------------------

Description

Normalizes intensity values to the median of each chip

Usage

```
normalizePerChipMedian(eset)
```

Arguments

eset	ExpressionSet containing chip intensity values to be normalized
------	---

Examples

```
#Normalize expression data in an Eset to the median
## Not run:
normedDataEset <- normalizePerChipMedian(eset)

## End(Not run)
```

outputAnnotatedDataMatrix

Outputs a tab delimited file from an ExpressionSet

Description

Takes an ExpressionSet and outputs a tab delimited file containing feature annotation to the left and hyb specific expression/flag data to the right

Usage

```
outputAnnotatedDataMatrix(eset, exptname, stage, dataElement)
```

Arguments

eset	ExpressionSet containing the matrix and annotation to output
exptname	a string containing the name of the experiment. Used to build file name
stage	a string containing the stage of the data in the matrix e.g. normalized
dataElement	a string containing the name of the data element in the ExpressionSet to be output

Examples

```
#Write out an annotated tab delimited file for the normalized data
## Not run:
outputAnnotatedDataMatrix(normedEset, "MyMicroArrayExpt", "Median_Normalized", "exprs")
## End(Not run)
```

panelCor

Add Pearson Correlation value to plots

Description

Adds a pearson correlation value to the scatter plots

Usage

```
panelCor(x,y, digits=2, prefix="r=")
```

Arguments

x	matrix of x values
y	matrix of y values to correlate with x
digits	number of digits to display
prefix	The string prefix that should be display on the scatterplot panel

Examples

```
#Calculate the median of a matrix omiting NAs
## Not run:
panelCor(x,y, digits=2, prefix="r=")

## End(Not run)
```

 parseRawData

Parse raw data files to create an ExpressionSet

Description

Loads all the gpr scanner output files in a particular directory and returns an ExpressionSet of the hybridizations in a MiChip experiment

Usage

```
parseRawData(datadir=".", pat="gpr")
```

Arguments

datadir	a directory containing one or my files of scanner output from MiChip hybridizations
pat	a string containing the three letter extension of the scanner output files

Examples

```
## Not run:
## Load all *.gpr files in current directory
parseRawData(datadir=".", pat="gpr")

## Load all *.gpr files in a specified directory , windows
parseRawData(datadir="c:\\mydata\\grpdata\\expt1\\", pat="gpr")

## Load all *.gpr files in a specified directory, linux
parseRawData(datadir="/home/myuser/gprdata/extp1/", pat="gpr")

## End(Not run)
```

plotIntensitiesScatter

Plot pairwise intensity scatter

Description

Creates a pairwise set of scatter plots from a data matrix and writes it out to file

Usage

```
plotIntensitiesScatter(dmat, controls=NULL, exptname, maintitle)
```

Arguments

dmat	matrix containing data from an experiment to be plotted
controls	matrix of row numbers containing control data to be plotted in a different colour
exptname	Name of the experiment, used for build the filename
maintitle	String used to build the maintitle of the graph

Examples

```
#Plot the pairwise intensities from myexpt
## Not run:
plotIntensitiesScatter(dmat, NULL, "MyExpt", "Median_Normalized")

## End(Not run)
```

removeUnwantedRows

Removes unwanted rows from data matrix

Description

Due to the requirements of spotting the chips, some of the spots are empty. Others contain controls or features from another species that may not be wanted in the analysis. This method removes them

Usage

```
removeUnwantedRows(rawData, filters)
```

Arguments

rawData	ExpressionSet containing matrix of data to be filtered
filters	list of strings to be filtered from annotation gene name column

Examples

```
#Removes empty and control spots from data matrix
## Not run:
filters=c("empty", "control")
filteredData <- removeUnwantedRows(rawData, filters)
## End(Not run)
```

```
returnAnnotatedDataMatrix
```

returns and annotated data matrix from an ExpressionSet

Description

Takes an ExpressionSet and returns a data matrix of feature annotation to the left and hyb specific expression/flag data to the right

Usage

```
returnAnnotatedDataMatrix(eset, dataElement)
```

Arguments

eset	ExpressionSet containing the matrix and annotation to output
dataElement	a string containing the name of the data element in the ExpressionSet to be output

Examples

```
#Write out an annotated tab delimited file for the normalized data
## Not run:
returnAnnotatedDataMatrix(normedEset,"exprs")
## End(Not run)
```

```
setIntensityCutoff
```

Sets a cutoff for the lowest intensity value

Description

Any value less than the cutoff value will be set to NA. This allows near background intensity values to be excluded

Usage

```
setIntensityCutoff(dmat, intensityCutoff)
```

Arguments

`dmat` matrix of intensity values to which the cutoff value is applied
`intensityCutoff` value of lowest acceptable intensity value in the experiment

Examples

```
#Set all the values under 50 in a matrix to NA
## Not run:
dmatOver50 <- setIntensityCutoff(dmat, 50)

## End(Not run)
```

<code>standardRemoveRows</code>	<i>Removes a standard list of features for MiChip processing</i>
---------------------------------	--

Description

Removes all empty spots, control spots, U6 RNA, non human spots from an ExpressionSet in the standard fashion. A wrapper for `removeUnwantedRows`

Usage

```
standardRemoveRows(rawData)
```

Arguments

`rawData` ExpressionSet containing the matrix to be filtered

Examples

```
#Filter standard rows from an ExpressionSet
## Not run:
myfilterdEset <-standardRemoveRows(rawData)

## End(Not run)
```

summarizeIntensitiesAsMedian

Summarizes the probe intensity as median of replicates spotted

Description

As the probes are spotted onto the in quaduplet or duplicate the values have to be combined in some way. This function takes the median of the intensities for the spots. Effectively the mean for duplicates. If less than half of the spots are present an NA is added

Usage

```
summarizeIntensitiesAsMedian(eset,minSumlength=0, madAdjust=FALSE)
```

Arguments

eset	ExpressionSet containing probe intensity data to be summarized
minSumlength	The lowest acceptable length of the matrix to calculate a median
madAdjust	if TRUE then summarized data will be filtered according to the MAD median absolute deviation and set to NA if the median is less than MAD

Examples

```
#Calculate the median of a matrix omiting NAs
## Not run:
summarizeIntensitiesAsMedian(eset,minSumlength=0,madAdjust=TRUE)

## End(Not run)
```

workedExampleMedianNormalize

Worked Example of MiChip Processing

Description

Loads a set of hybridizations into a matrix and them proceeds to filter, summarize and median normalize them

Usage

```
workedExampleMedianNormalize(exptname, intensityCutoff=0, datadir=".", minSumlength, madAdjust = FALSE)
```

Arguments

exptname	string indicating the name of the experiment
intensityCutoff	The intensity value for accepting the spots intensity value in the experiment
datadir	The directory where hybridization files are found.
minSumlength	Minimum exceptable number of values to summarize intensity value.
madAdjust	if TRUE then summarized data will be filtered according to the MAD median absolute deviation and set to NA if the median is less than MAD

Examples

```
#Normalize data in the current directory to the median per chip
datadir <- system.file("extdata", package="MiChip")
myNormedEset <-workedExampleMedianNormalize("MyExpt", intensityCutoff=0, datadir, minSumlength=0, madAdjust=TRUE)
```

```
workedExampleNotNormalizedData
```

Worked Example of MiChip Processing

Description

Loads a set of hybridizations into a matrix and then proceeds to filter and summarize these data

Usage

```
workedExampleNotNormalizedData(exptname, intensityCutoff=0, datadir=".", minSumlength, madAdjust = FALSE)
```

Arguments

exptname	string indicating the name of the experiment
intensityCutoff	The intensity value for accepting the spots intensity value in the experiment
datadir	The directory contain data from the experiment
minSumlength	Minimum exceptable number of values to summarize intensity value.
madAdjust	if TRUE then summarized data will be filtered according to the MAD median absolute deviation and set to NA if the median is less than MAD

Examples

```
#Summarizes the data in the current directory
## Not run:
mySummarizedEset <-workedExampleNotNormalizedData("MyExpt", intensityCutoff=0, datadir=".", minSumlength=0, madAdjust=TRUE)

## End(Not run)
```

Index

- * **aplot**
 - boxplotData, [2](#)
 - boxplotDataNoFile, [3](#)
 - panelCor, [6](#)
- * **datagen**
 - myForgivingMedian, [4](#)
 - naOmitMedian, [5](#)
 - summarizeIntensitiesAsMedian, [11](#)
- * **data**
 - correctForFlags, [3](#)
 - removeUnwantedRows, [8](#)
 - setIntensityCutoff, [9](#)
 - standardRemoveRows, [10](#)
- * **design**
 - workedExampleMedianNormalize, [11](#)
 - workedExampleNotNormalizedData, [12](#)
- * **documentation**
 - MiChip, [4](#)
- * **file**
 - outputAnnotatedDataMatrix, [6](#)
 - parseRawData, [7](#)
 - returnAnnotatedDataMatrix, [9](#)
- * **hplot**
 - plotIntensitiesScatter, [8](#)
- * **models**
 - normalizePerChipMedian, [5](#)

boxplotData, [2](#)
boxplotDataNoFile, [3](#)

correctForFlags, [3](#)

Introduction (MiChip), [4](#)

MiChip, [4](#)
myForgivingMedian, [4](#)

naOmitMedian, [5](#)
normalizePerChipMedian, [5](#)

outputAnnotatedDataMatrix, [6](#)

panelCor, [6](#)
parseRawData, [7](#)
plotIntensitiesScatter, [8](#)

removeUnwantedRows, [8](#)
returnAnnotatedDataMatrix, [9](#)

setIntensityCutoff, [9](#)
standardRemoveRows, [10](#)
summarizeIntensitiesAsMedian, [11](#)

workedExampleMedianNormalize, [11](#)
workedExampleNotNormalizedData, [12](#)