

# Package ‘FISHalyseR’

October 14, 2021

**Type** Package

**Title** FISHalyseR a package for automated FISH quantification

**Version** 1.26.0

**Date** 2015-04-08

**Author** Karesh Arunakirinathan <akaresh88@gmail.com>, Andreas Heindl  
<andreas.heindl@icr.ac.uk>

**Maintainer** Karesh Arunakirinathan <akaresh88@gmail.com>, Andreas  
Heindl <andreas.heindl@icr.ac.uk>

**Description** FISHalyseR provides functionality to process and analyse  
digital cell culture images, in particular to quantify FISH  
probes within nuclei. Furthermore, it extract the spatial  
location of each nucleus as well as each probe enabling spatial  
co-localisation analysis.

**VignetteBuilder** knitr

**License** Artistic-2.0

**Depends** EBImage,abind

**Suggests** knitr

**biocViews** CellBiology

**git\_url** <https://git.bioconductor.org/packages/FISHalyseR>

**git\_branch** RELEASE\_3\_13

**git\_last\_commit** 1ac9d07

**git\_last\_commit\_date** 2021-05-19

**Date/Publication** 2021-10-14

## R topics documented:

analyseParticles . . . . .	2
calculateMaxEntropy . . . . .	3
calculateThreshold . . . . .	4
computeIlluminationCorrection . . . . .	5
processFISH . . . . .	5

**analyseParticles**      *Analyse*

### Description

Cleans a given binary image according to area criteria specified by the user.

### Usage

```
analyseParticles(Image,MaxSize,MinSize, isMask)
```

### Arguments

Image	Binary image
MaxSize	Maximum area allowed for objects
MinSize	Minimum area allowed for objects
isMask	In case isMask=1, the function assumes that the binary images contains nuclei. Nuclei with an area smaller than MaxSize and greater than MinSize will be removed. If isMask=0, the function assumes that the binary images contains probes and subsequently probes with an area smaller than MinSize or larger than MaxSize are removed

### Value

Returns a labeled image

### Author(s)

Karesh Arunakirinathan

### Examples

```
f = system.file( "extdata", "SampleFISHgray.jpg", package="FISHalyseR")
img = readImage(f)

anaImg <- analyseParticles(img, 20000, 1000,0)
## anaImg contains now the cleaned-up image
```

---

calculateMaxEntropy    *Max Entropy thresholding*

---

## Description

The function converts a grayscale image to a binary image by computing a threshold using the Max Entropy method.

## Usage

```
calculateMaxEntropy(Image)
```

## Arguments

Image	grayscale image
-------	-----------------

## Details

Max Entropy thresholding can be used to detect the signals of probes in FISH cell culture images.

## Value

The function returns the threshold value

## Author(s)

Karesh Arunakirinathan

## References

J.N KANPUR, P.K SHAOO, A.K.C WONG: A New Method for Gray-Level picture thresholding Using the Entropy of the Histogram. In COMPUTER VISION, GRAPHICS AND IMAGE PROCESSING,1985 p 273-285

## See Also

calculateThreshold

## Examples

```
f = system.file( "exdata", "SampleFISHgray.jpg", package="FISHalyseR")
img = readImage(f)

t = calculateMaxEntropy(img)

## Threshold grayscale image using the value computed by the Max Entropy method
img[img<t] <- 0
img[img>=t] <- 1
```

---

<code>calculateThreshold</code>	<i>Compute threshold using Otsu's method</i>
---------------------------------	--

---

## Description

Computes the binary image of a grayscale image by using Otsu thresholding

## Usage

```
calculateThreshold(Image)
```

## Arguments

<code>Image</code>	grayscale image
--------------------	-----------------

## Details

The function computes a binary image using Otsu's method.

## Value

`calculateThreshold` returns the threshold value

## Author(s)

Karesh Arunakirinathan

## References

Nobuyuki Otsu: A threshold selection method from grey level histograms. In: IEEE Transactions on Systems, Man, and Cybernetics. New York 9.1979, S.62-66. ISSN 1083-4419

## See Also

`calculateMaxEntropy`

## Examples

```
f = system.file( "extdata", "SampleFISHgray.jpg", package="FISHalyseR")
img = readImage(f)

t = calculateThreshold(img)

##Threshold image using the value computed via Otsu's method
img[img<t] <- 0
img[img>=t] <- 1
```

---

computeIlluminationCorrection  
*Multidimensional Illumination Correction*

---

**Description**

Function to compute the multidimensional illumination correction (MDIC) using a stack of images

**Usage**

```
computeIlluminationCorrection(Images, pattern='*', AmountOfFile=6)
```

**Arguments**

Images	Directory containing the images
pattern	Filenames have to match the pattern specified here
AmountOfFile	Limit the amount of files used to compute the illumination gradient

**Value**

```
computeIlluminationCorrection  
return the image containing the illumination background
```

**Author(s)**

Andreas Heindl

**Examples**

```
illuCorrection = dirname(system.file( "extdata", "SampleFISHillu.jpg", package="FISHalyseR"))
```

---

processFISH                   *FISHalyseR - Automated fluorescence *in situ* hybridisation quantification in R*

---

**Description**

Function to automatically quantify FISH probes in cell-culture images.

**Usage**

```
processFISH(combinedImg, wriedir, bgCorrMethod = list(1, 100), channelSignals = NULL,  
            channelColours = NULL, sizeNucleus = c(5, 15000), sizeProbe = c(5, 100),  
            gaussigma = 20, outputImageFormat = ".png")
```

**Arguments**

combinedImg	Composite image of all available channels
writedir	Traget directory for output files
bgCorrMethod	Specifies the method used to correct for uneven background. Accepts only list types. Currently, four different methods are available: (1) Gaussian blurring, (2) Illumination correction image provided by the user, (3) multidimensional illumination correction (using a stack of images). In case no illumination correction should be applied, pass an empty list to the function
channelSignals	List of images containing the FISH probe
channelColours	List of colour vectors for each single channel
sizeNucleus	Minimum and maximum area (in pixel) of probes to be considred for further analysis
sizeProbe	Minimum and maximum area (in pixel) of probes to be considered for further analysis
gaussigma	Sigma of Gaussian used to blur the image
outputImageFormat	File format for the output image

**Value**

processFISH does not return any value

**Author(s)**

Karesh Arunakirinathan, Andreas Heindl

**See Also**

computeIlluminationCorrection, analyseParticles

**Examples**

```
## Specify illumination correction image
illuCorrection = system.file( "extdata", "SampleFISHillu.jpg", package="FISHalyseR")

## Composite image containing available channels
combinedImage <- system.file( "extdata", "SampleFISH.jpg", package="FISHalyseR")

## Single FISH channels containing the probe signals
red_Og   <- system.file( "extdata", "SampleFISH_R.jpg", package="FISHalyseR")
green_Gn <- system.file( "extdata", "SampleFISH_G.jpg", package="FISHalyseR")

## Output directory
writedir = paste(tempdir(),sep='')

## Use provided illumination correction image
bgCorrMethod = list(2,illuCorrection)
```

```
## Colour vector for three different probe channels (red, green and blue)
channelColours = list(R=c(255,0,0),G=c(0,255,0))

## Add probe channels to list
channelSignals = list(red_Og,green_Gn)

## Minimum and maximum area allowed for nuclei respectively probes
sizecell = c(1000,20000)
sizeprobe= c(5,20)

## Call processFISH with the specified parameters
processFISH(combinedImage,writedir,bgCorrMethod,channelSignals,
            channelColours,sizecell,sizeprobe)
```

# Index

```
* misc, bwlabel
    analyseParticles, 2
* misc
    calculateMaxEntropy, 3
    calculateThreshold, 4
    computeIlluminationCorrection, 5
    processFISH, 5

    analyseParticles, 2
    calculateMaxEntropy, 3
    calculateThreshold, 4
    computeIlluminationCorrection, 5

    processFISH, 5
```